

LAB325 Partial Thromboplastin Time (PTT)

BMC, CRH, LFH, OCH, FTT, FOX

Specimen Type

Platelet Poor Plasma

Preferred Container

Plastic blue top tube (3.2% sodium citrate)

Volume Required

2.7 mL

Storage Requirements and Transportation Needs

If transportation to the laboratory is to be delayed for more than 1 hour:

Prepare Platelet Poor Plasma:

Centrifuge tube at 1500 g for 10 minutes.

Use a plastic transfer pipette to remove plasma (staying away from the buffy coat layer) and transfer top 2/3 of plasma to a plastic aliquot tube.

Centrifuge this aliquot tube at 1500 g for 10 minutes.

Transfer top 3/4 of plasma (do not disturb button at bottom of tube) into another plastic aliquot tube.

Label this tube with patient information and a PLASMA sticker.

Freeze plasma immediately in a -15° to -25°C freezer until ready to transport.

Plasma must be transported to Laboratory frozen.

Causes for Rejection

Clotted specimens

Mislabeled or unlabeled specimens

Collection in any tube other than 3.2% plastic sodium citrate blue top tube

Less than 90% filled

Received greater than 4 hours old unless separated from cells and frozen at -15° to -25°C

Hemolysis > 500 mg/dL (>4+)

Whole blood or Plasma samples received refrigerated (2-8° C) or on ice

Reference Values

Population	Units	Reference Range	High Critical Value	Result Comment
Adult populations without heparin therapy	Seconds	27.1 – 39.1	≥ 69.2	Warfarin and Direct oral anticoagulants may elevate aPTT.
Suggested range for heparin therapy (0.3 – 0.7 U/mL)	Seconds	59.6 – 82.3	> 93.7	

It should be noted that children run 2-3 seconds above the adult PTT reference range without demonstrable coagulation disorder.

Available STAT

Yes

Methodology

Turbidimetric

CPT Code

85730