BMC

Performing Location(s)

BMC

Specimen Type

Serum

Preferred Container

Gold or Red top tube

Minimum Volume to Submit for Testing

2 mL Serum

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Does not need to be sent from outreach sites frozen. Can be sent spun and refrigerated in primary tube.

Sample tubes should be centrifuged within 2 hours of collection followed by transfer of the serum (red) to a labeled plastic, aliquot tube.

Gold tops do not need to be aliquoted

See Causes of Rejection for temperature requirements.

Transportation Needs

Deliver specimen to laboratory within 1 hour of collection if unspun.

If not possible see "Storage Requirements".

Causes for Rejection

Mislabeled or unlabeled specimen

Less than 50% draw for vacutainer tubes

- > 5 months refrigerated
- > 3 hours room temperature
- > 5 months frozen

Hemolysis > 500 mg/dL HgB

Lipemia > 1500 mg/dL

Trig Bilirubin > 20 mg/dL Bili

EDTA, Heparin, NaF, Citrate, Oxalate

TP < 4.2 or > 12.1 g/dL

Reference Values

Interpretation of Total and Free PSA Results

Probability (%) of Detecting Prostate Cancer: Men with Non-Suspicious DRE Results (Parentheses indicate Binomial 95% Confidence Intervals) Total PSA Reference Range: 0.1 - 4.0 ng/mL @Total PSA ... and %Free PSA of ... the %Probability of Detecting Prostate Cancer for Age Group 50-59 60-69 **70**+ 70.3 4.0-10.0 ng/mL < 10% 45.3 58.0 (33.8 - 57.3)(46.5 - 68.9)(53.0 - 84.1)11 - 19%22.5 30.3 38.0 (15.1 - 31.4)(23.6 - 37.7)(26.7 - 50.3)Greater than or 0.0 25.0 28.3 (N/A) equal to 20% (12.7 - 41.2)(16.0 - 43.5)

Reflex Testing

Prostate Cancer Prevalence (%):

If the Total Prostatic Antigen (PSA) is less than 4.0ng/mL, no further testing will be performed; only CPT 84153 will be billed.

30.3

37.2

42.9

If the Total Prostatic Antigen (PSA) is greater than 10.0ng/mL, no further testing will be performed; only CPT 84153 will be billed.

If the Total Prostatic Antigen (PSA) is greater than 4.0ng/mL and is less than 10.0ng/mL, a Free PSA will be performed and the %Free PSA will be calculated; CPT's 84153 and 84154 will be billed.

Available STAT

No

Methodology

Beckman AU/Dx Series

CPT Code

84153