LAB9918 BUN (NovaBio)

BMC, CRH, LFH, OCH, FTT, FOX



🔔 Important Note

LAB ORDER ONLY - BUN Prime Plus

Performing Location(s)

BMC, CRH, LFH, OCH, FTT, FOX

Specimen Type

Whole Blood

Preferred Container

Lithium Heparin Tube

Minimum Volume to Submit for Testing

2.0 mL Whole Blood

Transportation Needs

Do not centrifuge.

Deliver within 30 minutes after collection at RT.

Stable 24 hours at RT.

Causes for Rejection

Mislabeled or unlabeled specimen Less than 50% draw for vacutainer tubes

Clotted or Centrifuged

Specimen > 24 hours at room temperature

All other tubes or anticoagulant other than lithium heparin whole blood

Reference Values

Component	Reference Range
BUN	7-17 mg/dL

Available STAT

Yes

Methodology

NovaBio Medical Prime Plus

CPT Code

84520