BMC

Performing Location(s)

BMC

Specimen Type

Serum or Plasma

Preferred Container

Gold top tube

Alternate Container

Red top tube, Green top tube

Minimum Volume to Submit for Testing

 $3 \, mL$

Storage Requirements

Sample tubes should be centrifuged within 1 hour of collection followed by transfer of the serum (red) to a labeled plastic, aliquot tube.

Gold tops do not need to be aliquoted

Transportation Needs

Deliver specimen to laboratory within 1 hour of collection if unspun.

If not possible see "Storage Requirements".

Causes for Rejection

Mislabeled or unlabeled specimen

Less than 50% draw for vacutainer tubes

- > 48 hours refrigerated
- > 8 hours room temp
- > 6 months frozen

No Add-Ons

Performed only on fresh specimens

Reference Values

Population	Units	Reference Range
Female	mIU/mL	0-6
Male	mIU/mL	0-3

Available STAT

No

Methodology Immunoassay

CPT Code 84702