BMC, CRH, LFH, HRK, OCH, FTT, FOX, FCC

# **Performing Location(s)**

BMC, CRH, LFH, HRK, OCH, FTT, FOX, FCC

# **Specimen Type**

Serum or Plasma

### **Preferred Container**

Gold or Red top tube

### **Alternate Container**

Green top tube

## **Minimum Volume to Submit for Testing**

0.5 mL Serum or Plasma Multiple test can be performed from this volume. For question please contact the lab at 547-3975.

# **Storage Requirements**

Sample tubes should be centrifuged within 2 hours of collection Gold tops DO NOT ALIQUOT Aliquot Red tops after centrifuging See Causes of Rejection for temperature requirements.

### **Transportation Needs**

Deliver specimen to laboratory within 1 hour of collection if unspun. If not possible see "Storage Requirements".

## **Causes for Rejection**

Mislabeled or unlabeled specimen Less than 50% draw for vacutainer tubes > 7 days refrigerated > 7 days room temperature > 3 months frozen Hemolysis > 500 mg/dL HgB Lipemia > 1000 mg/dL Trig Bilirubin > 8 mg/dL Bili

### **Reference Values**

Population	Units	Reference Range
All Populations	mmol/L	98-110

# **Available STAT**

Yes

Methodology Beckman AU/Dx Series

# **CPT Code**

82435