

LAB712A Creatinine, 12 Hour Urine

BMC

Performing Location(s)

BMC

Specimen Type

12 hour urine

Preferred Container

24 Hour Urine Container with No Preservative
Keep refrigerated during collection

Minimum Volume to Submit for Testing

Entire collection
At least 10 mL urine needs to be submitted to testing lab.

Transportation Needs

Deliver specimen to the laboratory at the end of collection.
If transportation is expected to exceed 1 hour, provisions must be made to keep the urine chilled/

Storage Requirements

Keep urine refrigerated during and after the collection period.

Causes for Rejection

Mislabeled or unlabeled specimen
Greater than 10% of total volume lost
Inappropriate urine preservative used
> 2 days room temperature
> 6 days refrigerated
Specimen not refrigerated during and after collection

Reference Values

Population	Units	Reference Range
Female	mg/12 hour	0.3-0.8
Male	mg/12 hour	0.3-0.8

Available STAT

No

Methodology

Beckman AU/Dx Series

CPT Code

82570