

LAB384 Creatinine, Random Urine

BMC

Performing Location(s)

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see "Storage requirements"

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated

Causes for Rejection

Mislabeled or unlabeled specimen

Inappropriate urine preservative used

Specimen not refrigerated during or after collection

Hemolysis

> 2 days room temperature

> 6 days refrigerated

Reference Values

Population	Units	Reference Range
Male	mg/dL	40.0-370.0
Female	mg/dL	15.0-310.0

Additional Information

Note: If a urine protein to urine creatinine ration is desired, A Protein, Total: Urine, Random (LIS order code: RUP) must also ordered and a Urine Protein: Urine Creatinine Ratio will be calculated.

Available STAT

No

Methodology

Beckman AU/Dx Series

CPT Code

82570