

LAB21068 Dialysis, 24 Hour Urine Profile

BMC

Performing Location(s)

BMC

Specimen Type

24 Hour Urine

Preferred Container

24 Hr Urine Container with No Preservative

Alternate Container

No alternate container.

Minimum Volume to Submit for Testing

Entire collection

Storage Requirements

Keep urine refrigerated during the collection period.

Transportation Needs

Deliver specimen to the laboratory at the end of collection.

Causes for Rejection

Mislabeled or unlabeled specimen

Greater than 10% of total volume lost

Inappropriate urine preservative used

Specimen not refrigerated

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

81050, 84540, 82570