Fox

Specimen Type

Serum

Preferred: Serum Gel Acceptable: Red Top

Specimen Stability

Refrigerated

Reject Due To

Lipemic, hemolyzed or microbial contaminated samples may give poor results and should not be used.

Undiluted samples at room temperature > 8 hours.

Reference Values

Interpretation	Value
Negative	< 10
Equivocal	10 - 15
Positive	> 15
Reference Range: > 15 IU/mL	

Performing Lab

Fox

Methodology

Fluorescence Enzyme Immunoassay (FEIA)

CPT Code

86225