

LAB3191 ENA Screen (Symphony IgG) - No Reflex

Fox

Specimen Type

Serum

Preferred: Serum Gel

Acceptable: Red Top

Specimen Stability

Refrigerated

Reject Due To

Lipemic, hemolyzed or microbial contaminated samples may give poor results and should not be used.

Undiluted samples at room temp > 8 hours.

Reference Values

Interpretation	Value
ENA	
Negative	< 0.7
Equivocal	0.7 – 1.0
Positive	> 1.0
Reference Range: ≤ 1.0 Ratio	
dsDNA	
Negative	< 10
Equivocal	10 – 15
Positive	> 15
Reference Range ≤ 15 IU/mL	

Performing Lab

Fox

Methodology

Fluorescence Enzyme Immunoassay (FEIA)

CPT Codes

86038