BMC, FOX, HLS, FCC

Performing Location(s)

MIB, FOX, HLS, FCC

Specimen Type

Serum

Preferred Container

Gold or Red top tube

Alternate Container

Green top tube

Minimum Volume to Submit for Testing

1 ml Serum

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Transportation Needs

Deliver specimen to laboratory within 1 hour of collection if unspun.

If not possible see "Storage Requirements".

Storage Requirements

Sample tubes should be centrifuged within 2 hours of collection

Gold tops DO NOT ALIQUOT

Aliquot Red tops

See Causes of Rejection for temperature requirements.

Causes for Rejection

Mislabeled or unlabeled specimen

Less than 50% draw for vacutainer tubes

- > 7 days refrigerated
- > 24 hours room temperature
- > 12 months frozen

Hemolysis > 500 mg/dL HgB

Lipemia > 400 mg/dL Trig

Bilirubin > 40 mg/dL Bili

Creatinine > 30 mg/dL

Urea > 500 mg/dL

Uric acid > 20 mg/dL

Reference Values

Population	Units	Reference Range
Female	ng/mL	10-158
Male	ng/mL	16-243
Unknown Sex	ng/mL	10-243

Available STAT

Yes

Methodology Beckman AU/DX Series

CPT Code

82728