

LAB22222 Flow Cytometry, Other Source

URMC

Additional Test Codes

URMC: FLOWC

FLOWB Depending on source

Specimen Type

Body fluids, except CSF (For CSF, please order LAB22225).

Tissue - RPMI (Transport media tubes), sterile saline also acceptable. Specimens must be received within 48 hours of collection.

Special Requirements

Must have a **completed and signed** URMC requisition sent with specimens.

Place **specimens, requisition, CBCA results, and billing** in URMC shipping container.

Place bag in URMC pick-up area in CLP and enter information into the log sheet.

If specimens are collected at sites other than MIB campus they MUST arrive by courier by 6pm. Courier arrives for pick up at 7pm.

Storage/Transport Temperature

Tissue: 2-8°C

Body Fluids- Room Temp

Causes of Rejection

Body fluids and tissue in formalin/fixed

Frozen or warmer than 37°C

CPT Code

88184 (for 1 marker)

88187 (2-8 markers);

88188 (9-15 markers)

88189 (Greater than 16 markers 88189).

Number of markers tested is dependent on the diagnosis.