

LAB9917 Glucose (NovaBio)

BMC, CRH, LFH, OCH, FTT, FOX

**Important Note**

LAB ORDER ONLY - GLUCOSE Prime Plus

Performing Location(s)

BMC, CRH, LFH, OCH, FTT, FOX

Specimen Type

Whole Blood

Preferred Container

Lithium Heparin Tube

Minimum Volume to Submit for Testing

2.0 mL Whole Blood

Transportation Needs

Do not centrifuge.

Deliver within 30 minutes after collection at room temperature.

Stable 1 hour at room temperature.

Causes for Rejection

Mislabeled or unlabeled specimen

Less than 50% draw for vacutainer tubes

Clotted or Centrifuged

Specimen > 1 hour at room temperature.

All other tubes or anticoagulant other than lithium heparin whole blood

Reference Values

Component	Reference Range	Low Critical	High Critical
Glucose, Random 0 to 7 Days	45-100 mg/dL	<45 mg/dL	>450 mg/dL
Glucose, Random >7 Days	70-139 mg/dL	<40 mg/dL	>450 mg/dL

Available STAT

Yes

Methodology

NovaBio Medical Prime Plus

CPT Code

82947