

LAB186 Glucose, Body Fluid

BMC

Performing Location(s)

BMC

Specimen Type

Body Fluid

Preferred Container

Red top tube

Alternate Container

Green top tube

Minimum Volume

0.5 mL Body Fluid

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to laboratory within 1 hour should be refrigerated.

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

Specimen too viscous to analyze

Clotted

Collected in wrong tube

> 2 hours at room temperature

> 4 days refrigerated

Reference Values

Fluid Type	Units	Reference Range
Synovial fluid	mg/dL	Not Established
Other than synovial fluid	mg/dL	Not Established

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

82945