

LAB21029 Glucose, Random Urine

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume

1 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

Specimen collected with HCl as preservative

Reference Values

Population	Units	Reference Range
All Populations	mg/dL	1 - 15

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

82945