BMC

Performing Location(s)

BMC

Specimen Type

Body Fluid

Preferred Container

Sterile screw top container

Alternate Container

Gray top tube

Minimum Volume to Submit for Testing

0.5 mL Body Fluid

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to laboratory within 1 hour should be refrigerated

Transportation Needs

Deliver specimen to the laboratory **ON ICE** within 1 hour of collection.

If possible see "Storage Requirements"

Causes for Rejection

Mislabeled or unlabeled specimen Not on ice within 1 hour of collection Specimen too viscous to analyze > 3 days refrigerated

Reference Values

Population	Units	Reference Range
All populations	U/L	No Established

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code