BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume to Submit for Testing

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated.

Transportation Needs

Deliver specimens to the laboratory within 1 hour of collection.

If not possible, see "Storage Requirements"

Causes for Rejection

Mislabeled or unlabeled specimen

Specimen collected with unacceptable preservative

- > 3 days room temperature
- > 3 days refrigerated

Reference Values

| Population | Units | Reference Range |
|-----------------|-------|-----------------|
| All Populations | mg/dL | 1.0 - 13.0 |

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

83735