

LAB694 Occult Blood, Stool (Single Card)

BMC, CRH, LFH, OCH, FTT, FOX

Specimen Type

Random Stool

Preferred Container

Hemoccult Sensa

Alternate Container

Stool Container

Volume Required

Single Card

Transportation Needs

Cards must be returned to the laboratory for testing within 21 days after preparation.

*Transport at room temperature

Storage Requirements

Store the card at room temperature

Causes for Rejection

Specimen is applied to wrong side of the card

Stool card received beyond expiration date printed on card

Stool card no received within 21 days of collection

Stool samples > 48 hours old if collected and stored in a stool container

Reference Values

| Population | Units | Reference Range |
|-----------------|-------|-----------------|
| All Populations | none | Negative |

Available STAT

Yes

CPT Code

82272