

LAB695 Occult Blood, Stool (Triple Cards)

BMC, CRH, LFH, OCH, FTT, FOX



Important Note

Only Orderable at: MIB GI LAB, MIB Medical Clinic, OSS Specialty, Cobleskill Specialty and Herkimer Specialty

Used to assist with ruling out upper GI Bleeds

Specimen Type

3 consecutive random stool samples within 21 days.

Preferred Container

Hemoccult Sensa

Alternate Container

Stool Container

Volume Required

Three Cards

Transportation Needs

Cards must be returned to the laboratory for testing within 21 days after preparation.

*Transport at room temperature

Storage Requirements

Store the card at room temperature

Causes for Rejection

Specimen is applied to wrong side of the card

Stool card received beyond expiration date printed on card

Stool card no received within 21 days of collection

Stool samples > 48 hours old if collected and stored in a stool container

Reference Values

Population	Units	Reference Range
All Populations	none	Negative

Available STAT

Yes

CPT Code

82270