

POC189 POC Glucose

Bassett Healthcare Network POC Testing

Performing Location(s)

Bassett Healthcare Network POC Testing. Only available at POC Approved Testing Locations.

Specimen Type

Capillary

Preferred Container

Capillary (non-neonate fingerstick or neonate heel stick)

Volume Required

1.2 uL

Additional Information

Fresh venous and arterial samples containing lithium heparin are acceptable samples if well mixed and tested immediately.

DO NOT use EDTA or sodium heparin.

Repeating Samples: Repeat test if results are not clinically expected.

Reference Values

Age	Normal	High Critical	Low Critical
<1 Week	45-100 mg/dL	>450 mg/dL	<45 mg/dL
≥1 Week	70-139 mg/dL	>450 mg/dL	<40 mg/dL

Confirmatory Testing

Venous Confirmation: A Laboratory confirmation is required for results < 10 (LO) or > 600 (HI).

Available STAT

Only available at POC Approved Testing Locations.

Methodology

Nova Stat Strip

CPT Code

82948