

## POC149 POC Hemoglobin

Bassett Healthcare Network POC Testing

### Performing Location(s)

Bassett Healthcare Network POC Testing. Only available at POC Approved Testing Locations.

### Specimen Type

Capillary

### Preferred Container

Wipe away first two or three drops.

### Volume Required

10 uL into HemoCue Microcuvette

Avoid bubbles and under-filling

### Reference Values

Population	Value
< 1 month	13.0 – 19.0 g/dL
1 – 5 months	11.0 – 17.0 g/dL
5 months – 2 years	11.0 – 14.0 g/dL
2 – 5 years	11.0 – 13.0 g/dL
6 – 12 years	11.0 – 15.0 g/dL
> 12 years male	11.5 – 18.0 g/dL
> 12 years female	11.5 – 15.5 g/dL
<b>Critical Value: &lt; 7.0 g/dL</b>	

### Confirmatory Testing

Venous Confirmation: A Laboratory confirmation is required for results <6.0.

### Additional Information

Collect samples from the middle finger ONLY.

Repeating Samples: All results < 7.0 must be repeated with a new fingerstick or confirmatory tube sent to the laboratory.

### Available STAT

Only available at POC Approved Testing Locations.

### Methodology

iSTAT