

POC550 POC Urinalysis (Clinitek)

Bassett Healthcare Network POC Testing

Performing Location(s)

Bassett Healthcare Network POC Testing. Only available at POC Approved Testing Locations.

Specimen Type

Urine

Preferred Container

Sterile container (required for culture)

Volume Required

5 mL

Macroscopic			
Parameter	Population	Units	Reference Range
Color	All populations	None	Yellow, straw or colorless
Appearance	All populations	None	Clear
Specific Gravity	All populations	None	1.003-1.024
pH	All populations	None	5.0-8.0
Protein	All populations	mg/dL	Negative
Glucose	All populations	mg/dL	Negative
Ketone	All populations	mg/dL	Negative
Bilirubin	All populations	None	Negative
Blood	All populations	None	Negative
Urobilinogen	All populations	Eu/dL	≤ 1.0
Nitrite	All populations	None	Negative
Leuk Esterase	All populations	None	Negative

Confirmatory Testing

Positive Bilirubin should be confirmed by laboratory testing.

Positive nitrites and/or positive leukocyte esterase in screening for UTI should follow with urine culture (LAB239)

Available STAT

Only available at POC Approved Testing Locations.

Methodology

Clinitek Status+

CPT Code

81003