

POC119 POC Urine, Microscopic

Bassett Healthcare Network POC Testing

Performing Locations

Bassett Healthcare Network POC Testing. Only available at POC Approved Testing Locations.

Specimen Type

Random urine

Preferred Container

Random Urine Container

Place one drop of specimen onto glass slide and apply cover slip

Minimum Volume for Testing

10 mL Urine

Reference Values

Microscopic			
Parameter	Population	Units	Reference Range
RBC's	All populations	#/hpf	0-2
WBC's	All populations	#/hpf	<6
Squamous Epithelium	All populations	#/hpf	<6
Transitional Epithelium	All populations	#/hpf	0-3
Bacteria	All populations	None	None
Crystals	All populations	#/hpf	None
Hyaline Casts	All populations	#/lpf	0-2
All Other Casts	All populations	#/lpf	None

Additional Information

Interpretation of formed elements can be subjective. Submit urine to the Clinical Laboratory for confirmation of cellular elements and special staining for fats or oval fat bodies.

Urine may be stored at room temperature for 1 hour.

Available STAT

Only available at POC Approved Testing Locations

Methodology

Microscope

CPT Code

81015