

Varies Peritoneal Dialysates Effluent Culture

BMC

Performing Location(s)

BMC

Specimen Type

Dialysate

Preferred Container

Sterile Container

Volume Required

5 mL

Storage Requirements

Incubate 35°C

Transportation Needs

Room Temperature

Reference Values

Population	Units	Reference Range
All Populations	(none)	No Growth

Includes

Gram stain

Aerobic Culture

Anaerobic Culture

Methodology

Aerobic and Anaerobic Culture

CPT Code

87070

87075

87205