

## **LAB434 Potassium, Random Urine**

---

BMC

### **Performing Location(s)**

BMC

### **Specimen Type**

Random Urine

### **Preferred Container**

Random Urine Tube

### **Minimum Volume to Submit for Testing**

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

### **Storage Requirements**

Specimens not delivered to laboratory within 1 hour should be refrigerated

### **Transportation Needs**

Deliver specimens to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

### **Causes for Rejection**

Mislabeled or unlabeled specimen

> 12 hours room temperature

> 45 days refrigerated

Collected with a preservative

### **Reference Values**

| <b>Population</b> | <b>Units</b> | <b>Reference Range</b> |
|-------------------|--------------|------------------------|
| All Populations   | mmol/L       | 12.0-62.0              |

### **Available STAT**

Yes

### **Methodology**

Beckman AU/Dx Series

### **CPT Code**

84133