**BMC** 

## **Performing Location(s)**

**BMC** 

# **Specimen Type**

Random Urine

#### **Preferred Container**

Random Urine Tube

## **Minimum Volume to Submit for Testing**

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

## **Storage Requirements**

Keep urine refrigerated

### **Transportation Needs**

Deliver specimen to the laboratory within 1 hour of collection.

If transportation is expected to exceed 1 hour, provisions must be made to keep the urine chilled.

## **Causes for Rejection**

Mislabeled or unlabeled specimen Bloody specimens Hemolysis Collected with preservative > 12 hours room temperature

> 2 days refrigerated

### **Reference Values**

Population	Units	Reference Range
All populations	mg/dL	<15

#### **Available STAT**

Yes

#### Methodology

Beckman AU/Dx Series

#### **CPT Code**

84156