## LAB21006 pH, Gastric

**BMC** 

# **Performing Location(s)**

BMC, FTT

# **Specimen Type**

Gastric Contents, vomitus or nasogastric

#### **Preferred Container**

Sterile container

#### **Alternate Container**

Random Urine Tube

## **Volume Required**

1 mL

## **Storage Requirements**

Testing must be completed within 2 hours of collection.

## **Transportation Needs**

Deliver to laboratory immediately.

Testing needs to be completed within 2 hrs of collection.

A STAT Courier will need to be arranged to have specimen delivered to the testing location.

## **Causes for Rejection**

Mislabeled or unlabeled specimens Insufficient quantity Grossly bloody

- > 2 hours at room temperature
- > 2 hours refrigerated

#### **Reference Values**

| Population      | Units | Reference Range |
|-----------------|-------|-----------------|
| All populations | None  | 1.0 - 3.0       |

## **Available STAT**

Yes

# Methodology

pH paper – hydrogen ion concentration

# **CPT Code**