

## LAB3197 RNP Antibodies, IgG

---

Fox

### Specimen Type

Serum

Preferred: Serum Gel

Acceptable: Red Top

### Specimen Stability

Refrigerated

### Reject Due To

Lipemic, hemolyzed or microbial contaminated samples may give poor results and should not be used.

Undiluted samples at room temp > 8 hours.

### Reference Values

Component	Negative	Equivocal	Positive
RNP70	< 7	7 – 10	> 10
U1RNP	< 5	5 – 10	> 10

### Performing Location

Fox

### Methodology

Fluorescence Enzyme Immunoassay (FEIA)

### CPT Codes

86235 × 2