

**LAB22190** Separate, Freeze and Hold (Plasma)

---

BMC

**Performing Location(s)**

BMC

**Specimen Type**

Plasma

**Preferred Container**

Green, Lavender, or Blue top tube

**Storage Requirements**

Freeze

**Transportation Needs**

Deliver specimen to laboratory within 1 hour of collection if unspun.

If not possible see "Storage Requirements".

**Causes for Rejection**

Mislabeled or unlabeled specimen

Less than 50% draw for vacutainer tubes

**CPT Code**

No charge associated with this test