

**LAB22145** Separate, Freeze and Hold (Serum)

---

BMC

**Performing Location(s)**

BMC

**Specimen Type**

Serum

**Preferred Container**

Gold or Red top Tube

**Storage Requirements**

Freeze

**Transportation Needs**

Deliver specimen to laboratory within 1 hour of collection if unspun.

If not possible see "Storage Requirements".

**Causes for Rejection**

Mislabeled or unlabeled specimen

Less than 50% draw for vacutainer tubes

**CPT Code**

No charge associated with this test