

LAB3196 Smith (Sm) Antibodies, IgG

Fox

Specimen Type

Serum

Preferred: Red top

Acceptable: Serum Gel

Transport Temperature

Refrigerated

Performing Location

Fox

Reject Due To

Lipemic, hemolyzed or microbial contaminated samples may give poor results and should not be used.

Undiluted samples at room temp > 8 hours.

Reference Values

Interpretation	Value
Negative	< 7
Equivocal	7 – 10
Positive	> 10

Methodology

Fluorescence Enzyme Immunoassay (FEIA)

CPT Codes

86235