

LAB3183 Thyroid Autoantibodies Profile

Fox

Specimen Type

Serum

Preferred: Serum Gel

Acceptable: Red Top

Specimen Stability

Refrigerated

Reject Due To

Lipemic, hemolyzed or microbial contaminated samples may give poor results and should not be used.

Undiluted samples at room temp > 8 hours.

Reference Values

Test	Negative	Equivocal	Positive
Thyroglobulin Antibody	< 40	40 – 60	> 60
Thyroperoxidase Antibodies	< 25	25 – 35	> 35

Performing Location

Fox

Methodology

Fluorescence Enzyme Immunoassay (FEIA)

CPT Codes

86800

86376