

LAB21003A Urea Nitrogen, Dialysate Fluid

BMC

Performing Location(s)

BMC

Specimen Type

Dialysate Fluid

Preferred Container

Sterile Container

Minimum Volume to Submit for Testing

0.5 mL Dialysate Fluid

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to laboratory within 1 hour should be refrigerated

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

Causes for Rejection

Mislabeled or unlabeled specimen

Specimen too viscous to analyze or clotted

Specimen frozen

Collected in wrong tube

> 2 hours at room temperature

> 4 days refrigerated

Reference Values

Population	Units	Reference Range
All Populations	mg/dL	Not established

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

84520