BMC

Performing Location(s)

BMC

Specimen Type

44 Hour Urine

Preferred Container

24 Hr Urine Container with No Preservative.

Alternate Container

No alternate container.

Minimum Volume to Submit for Testing

1 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Keep urine refrigerated during the collection period

Transportation Needs

Deliver specimen to the laboratory at the end of collection keeping specimen refrigerated.

Causes for Rejection

Mislabeled or unlabeled specimen Greater than 10% of total volume lost Inappropriate urine preservative used Specimen not refrigerated during and after collection

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

81050, 84540, 82570