BMC

Performing Location(s)

BMC

Specimen Type

Site dependent

Preferred Container

eSwab (see also Specific Specimen Requirements)

Alternate Container

Swab specimens will be accepted, but are suboptimal and should be submitted only if no alternative exists

Transportation Needs

Specimens for anaerobic cultures must be transported anaerobically to the laboratory immediately

Storage Requirements

Keep specimen at room temperature

Causes for Rejection

Mislabeled or unlabeled specimen Sample > 72 hours old

Reference Values

Population	Units	Reference Range
All populations	(none)	No Growth

Reflex Testing

If a portion of specimen is submitted in an anaerobic transport Vial or tube and, if the source is appropriate, an anaerobic culture will also be performed and billed (CPT: 87075).

Organism identification and susceptibility studies are performed as appropriate, according to established protocols that have been approved by the Clinical Laboratory Director. Call Microbiology, 547-3707, for details.

There may be additional charges for these identification and susceptibility procedures

Includes

An order for a wound culture normally includes:

Aerobic Culture (CPT: 87070)

Gram Stain if adequate specimen submitted (2 culturette, eSwab or specimen in sterile container) (CPT: 87205)

Available STAT

No

Methodology Aerobic (& Anaerobic) culture

CPT Code

87070, 87205, 87075 (see INCLUDES)