

LAB21110 Calcium-Phosphorus Product

BMC

Performing Location(s)

BMC

Specimen Type

Serum or Plasma

Preferred Container

Gold or Red top

Alternate Container

Green top

Transportation Needs

Deliver specimen to laboratory within 1 hour of collection if unspun.

If not possible see "storage Requirements".

Causes for Rejection

Mislabeled or unlabeled original specimen Less than 50% draw for Vacutainer tubes

Reference Values

Population	Units	Reference Range
All Populations	(none)	0-55 mg ² /dL ²

Additional Information

NOTE: This test CANNOT be ordered as a single test. A Calcium (CALC) and a Phosphorus (PHOS) must be included with all requests for a Calcium/Phosphorus Product (CAPHOS)

Additional Information:

A Calcium (CPT: 82330) and a Phosphorus (CPT: 84100) will be performed and billed. There is no charge for the CAPHOS calculation

Available STAT

Yes

Methodology

Calculated from Calcium (CALC) and a Phosphorus (PHOS)

CPT Code

See Additional Information

