

LAB374 Chloride, Random Urine

BMC

Performing Location(s)

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated.

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “storage requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

> 12 hours room temperature

> 24 hours refrigerated

Specimen collected with HCl as preservative

Reference Values

Population	Units	Reference Range
All Populations	mmol/L	55 - 125

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

82436