## LAB39031 Low Dose Opioid, Urine

### **URMC**

# **Specimen Type**

Random Urine

# **Volume Required**

5 mL urine with no additives or preservatives

### **Minimum Volume**

3 mL urine

### **Transportation Needs**

Refrigerated

## **Storage Requirements**

Ambient: 1 week Refrigerated: 1 week

Frozen for long term storage.

# **Additional Information**

Includes parent drug and/or metabolites. Parent drug with cut-offs:

Codeine (50 ng/mL)

Morphine (30 ng/mL)

Hydrocodone (10 ng/mL)

Hydromorphone (30 ng/mL)

Oxycodone (50 ng/mL)

Oxymorphone (20 ng/mL)

Fentanyl (10 ng/mL)

Buprenorphine (3 ng/mL)

Maps to URMC: COPS (43760), FENT (19877), UBUP (30425)

**URMC-Drugs of Abuse Confirmation/Qualitative - Opiates - Urine (COPS)** 

**URMC-Drugs of Abuse Screen - Fentanyl & Metabolite - Urine (FENT)** 

<u>URMC-Buprenorphine & Metabolites - Quantitation - Urine (UBUP)</u>

### Methodology

Liquid Chromatography-Tandem Mass Spectrometry

### **CPT Code**

G6058, G0431, G6056