BMC, CRH, LFH, HRK, OCH, FTT, FOX, FCC

### **Performing Location(s)**

BMC, CRH, LFH, HRK, OCH, FTT, FOX, FCC

# **Specimen Type**

Whole Blood

#### **Preferred Container**

Lavender Top

# **Minimum Volume to Submit for Testing**

3.0 ml whole blood

#### **Transportation Needs**

Deliver to lab within 12 hours at room temperature.

## **Storage Requirements**

Refrigerate up to 24 hours.

### **Causes for Rejection**

Mislabeled or unlabeled specimen Less than 50% draw for vacutainer tubes Clotted

#### **Reference Values**

Population	Reference Range (cells/uL)
Newborn - 1 month	1600 - 6115
1 month - 5 months	750 - 3452
5 months - 2 years	650 - 3255
2 years - 6 years	750 - 3452
6 years - 12 years	1600 - 5326
> 12 years	1500 - 7400
Note: The absolute N	eutrophil Count includes mature
neutrophil	s & band neutrophils

# **Reflex Testing**

When an absolute Neutrophil Count is ordered on blood and a WBC and Differential (automated or manual) has not been ordered, a WBC and Differential (CPT's 85048 and 85004) will be performed and billed.

When an absolute Neutrophil Count is ordered on a blood that already has a WBC and Differential (automated or manual) ordered, it will be added on at no extra charge.

# **Available STAT**

Yes

**Methodology** Flow cytometry

**CPT Code** 85048