**BMC** 

## **Performing Location(s)**

**BMC** 

## **Specimen Type**

Random Stool

#### **Preferred Container**

Stool sample in modified Carey Blair (ETM)



### **Required Volume**

1g or 5mL Stool

#### **Storage Requirements**

Ambient temperature- preferred for up to 4 days (15 - 25°C)

### **Transportation Needs**

Deliver to the laboratory as soon as possible

## **Causes for Rejection**

Mislabeled or unlabeled specimen
Stool in Proto-Fix, formalin, SAF, MIF, PVA
Unpreserved specimen > 2 hours at room temperature
Swab - Rectal, Stool
Frozen specimen
Improper collection
Bloody stool

Mucus in stool Non-stool sample

# **Reference Values**

Population	Units	Reference Range
All Populations	(none)	Negative

# Available STAT

No

# **Methodology** BioFire

# **CPT Code**

87328