

POC7 POC Urine Pregnancy

Bassett Healthcare Network POC Testing

Performing Location(s)

Bassett Healthcare Network POC Testing. Only available at POC Approved Testing Locations.

Specimen Type

Random urine

Preferred Container

Urine collection cup

Minimum Volume for Testing

5 mL Urine

Reference Values

Reference Values: Negative

Additional Information

The first morning specimen of urine usually contains the highest concentration of hCG and is the best sample.

If urine appears cloudy or turbid, centrifugation prior to testing is recommended.

Confirmatory Testing

Serum hCG can be ordered

Available STAT

Only available at POC Approved Testing Locations.

Methodology

Rapid monoclonal and polyclonal membrane immunoassay

CPT Code

81025