

POC112 POC Vaginal pH

Bassett Healthcare Network POC Testing

Performing Location(s)

Bassett Healthcare Network POC Testing. Only available at POC Approved Testing Locations.

Specimen Type

Vaginal secretion

Preferred Container

Nitrazine paper

Collection Method

Apply drop of specimen onto pH paper and read immediately.

Reference Values

Value	Interpretation
< 4.5	Normal
Membrane Rupture in Pregnant Women	
≥ 6.5	Consistent with ruptured membrane
≤ 6.0	Consistent with intact membrane bacterial vaginosis
< 4.5	Virtually excludes bacterial vaginosis
>4.5	Consistent with bacterial vaginosis

Available STAT

Only available at POC Approved Testing Locations

Methodology

Nitrazine paper

CPT Code

83986