

LAB427 Phosphorus, Random Urine

BMC

Performing Location(s)

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine tube

Volume Required

2 mL

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated.

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

Specimen more than 8 hours old not refrigerated

Reference Values

Population	Units	Reference Range
All Populations	mg/dL	20 - 60

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

84105