

LAB444 Sodium, Random Urine

BMC

Performing Location(s)

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume to Submit for Testing

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated

Transportation Needs

Deliver specimens to the laboratory within 1 hour of collection.

If not possible, see “Storage Requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

> 12 hours room temperature

> 45 days refrigerated

Collected with a preservative.

Reference Values

Population	Units	Reference Range
All Populations	mmol/L	20-110

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

84300