

## LAB200 Triglycerides, Body Fluid

---

BMC

### Performing Location(s)

BMC

### Specimen Type

Body Fluid

### Preferred Container

Red top tube

### Alternate Container

Green top tube

### Minimum Volume to Submit for Testing

0.5 mL Body Fluid

### Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

### Storage Requirements

Specimens not delivered to laboratory within 1 hour should be refrigerated

### Causes for Rejection

Mislabeled or unlabeled specimen

Specimen too viscous for analysis

Clotted

Oxalate or Sodium Fluoride

> 2 hours at room temperature

> 4 days refrigerated

### Reference Values

Population	Units	Reference Range
All Populations	gm/dL	Not Established

### Available STAT

Yes

### Methodology

Beckman AU/Dx Series

### CPT Code

82042