

LAB748 Urea Nitrogen, Random Urine

BMC

Performing Location(s)

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume to Submit for Testing

5 mL Urine

Multiple test can be performed from this volume.

For question please contact the lab at 547-3975.

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

> 12 hours room temperature

> 24 hours refrigerated

Collected with a preservative

Reference Values

Population	Units	Reference Range
All Populations	mg/dL	350 - 1000

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

84540