

LAB21038 Uric Acid, Random Urine

BMC

Performing Location(s)

BMC

Specimen Type

Random Urine

Preferred Container

Random Urine Tube

Minimum Volume to Submit for Testing

5 mL urine

Storage Requirements

Specimens not delivered to the laboratory within 1 hour should be refrigerated.

Transportation Needs

Deliver specimen to the laboratory within 1 hour of collection.

If not possible see “Storage Requirements”

Causes for Rejection

Mislabeled or unlabeled specimen

> 12 hours room temperature

>24 hours refrigerated

Greater than 10% volume lost

Collected with a preservative

Reference Values

Population	Units	Reference Range
All Populations	mg/dL	7.5 - 49.5

Available STAT

Yes

Methodology

Beckman AU/Dx Series

CPT Code

84560